

Carebidet Quick Guide 02

Getting Started

Preparing Carebidet for Use

Use this with the training video

This quick guide summarizes the procedure shown in the approved Carebidet training video. For complete operating instructions and cautions, use the Carebidet User Manual and the full video.



Quick steps

1. Open the Carebidet main-body cover.
2. Fill the clean-water tank with water and place it in the main body in the correct orientation.
3. Insert the Diaper Cup main-body connector into the side plate and push it in until the connection clicks securely into place.
4. Push the main-body connector into the tube connector until the clicking sound confirms the connection.
5. Plug the Diaper Cup power/electrical cable into the cable port next to the control panel.
6. Plug the Carebidet power cord into an electrical outlet and turn the unit on.
7. Confirm the display is operating normally and follow the on-screen prompt to select the required operation mode.

Important

Stop and re-check the setup if the individual has pain, skin pinching, leakage, an unresolved alert, or an insecure connection.