


Carebidet Quick Guide 03

Fitting & Daily Use

Wearing the Diaper Cup

Use this with the training video

This quick guide summarizes the procedure shown in the approved Carebidet training video. For complete operating instructions and cautions, use the Carebidet User Manual and the full video.



Quick steps

1. Prepare the Motion Control Wrap by opening the abdominal and waist bands.
2. Position the Diaper Cup within the wrap and secure the lower zipper so the lower portion of the cup is covered.
3. Assist the individual into the required position and bring the cup into place at the perineal area.
4. Center the private area within the Diaper Cup. Alternate left and right as needed so the body is positioned close to the center of the cup.
5. Pull up and secure the abdominal band to the waist band.
6. Check the leakage-prevention bands and the Diaper Cup position; make sure the cup is seated without pinching or folding the skin.
7. Confirm the tube connector is routed correctly and check the Carebidet control-panel status before use.

Important

Stop and re-check the setup if the individual has pain, skin pinching, leakage, an unresolved alert, or an insecure connection.