

Carebidet Quick Guide 04

Fitting & Daily Use

Changing Position

Use this with the training video

This quick guide summarizes the procedure shown in the approved Carebidet training video. For complete operating instructions and cautions, use the Carebidet User Manual and the full video.



Quick steps

1. Use the postural-change cushion when moving the individual between supine and lateral positions.
2. Place the cushion in the correct position before turning, then support the individual with an additional cushion if needed for comfort.
3. Keep the Diaper Cup and the supporting cushion level. Do not allow the cup to hang off the support or become only partially supported.
4. When turning to the side, maintain the Diaper Cup position and check that the body remains centered in the cup.
5. When returning toward supine, hold the lower portion of the Diaper Cup and open the upper zipper as shown in the training video.
6. After repositioning, re-check the cup, wrap, tubing, and skin for secure placement before resuming use.

Important

Stop and re-check the setup if the individual has pain, skin pinching, leakage, an unresolved alert, or an insecure connection.