

Carebidet Quick Guide 06

Fitting & Daily Use

Applying the Motion Control Wrap - Male

Use this with the training video

This quick guide summarizes the procedure shown in the approved Carebidet training video. For complete operating instructions and cautions, use the Carebidet User Manual and the full video.



Quick steps

1. Lay out the Motion Control Wrap and prepare the Diaper Cup and male module before positioning the individual.
2. Attach the male module to the Diaper Cup/Wrap assembly as shown in the training video.
3. Position the wrap under the individual and bring the Diaper Cup into place.
4. Turn the male module upside down and gently place the male genitalia into the module.
5. Rotate the male module back into the operating position and connect it to the Diaper Cup.
6. Secure the Motion Control Wrap around the individual and check that the module and cup are centered without pinching the skin.
7. Verify the tubing and connections before starting Carebidet.

Important

Stop and re-check the setup if the individual has pain, skin pinching, leakage, an unresolved alert, or an insecure connection.