

Carebidet Quick Guide 14

Cleaning & Routine Maintenance

Tube Disinfection

Use this with the training video

This quick guide summarizes the procedure shown in the approved Carebidet training video. For complete operating instructions and cautions, use the Carebidet User Manual and the full video.



Quick steps

1. Attach the Diaper Cup holder to the side of the main body and lock the lower part of the cup into the holder.
2. Replace the detachable cushion with the plastic cleaning cover.
3. Use the control panel to enter the Tube Cleaning menu, then change the mode to No. 2 Tube Disinfection.
4. Put the appropriate amount of the approved disinfectant into the disinfectant tube kit.
5. Press OK to start the Tube Disinfection cycle and allow it to complete.
6. When complete, remove the disinfectant tube kit from the plastic cleaning cover.
7. Disconnect the tube as shown, clean the sewage and water tanks, reinstall them, and close the main body.

Important

Stop and re-check the setup if the individual has pain, skin pinching, leakage, an unresolved alert, or an insecure connection.