

Carebidet Quick Guide 15

Storage

Long-Term Storage

Use this with the training video

This quick guide summarizes the procedure shown in the approved Carebidet training video. For complete operating instructions and cautions, use the Carebidet User Manual and the full video.



Quick steps

1. Pull the sewage-tank detaching button back and separate the tube connection.
2. Lift the main-body connector to separate the tube from the main body.
3. Connect the drain hose to the water outlet in the main body and place the drain hose into the water tank.
4. On the control panel, select Cleaning Settings > 1-015 Manual Drainage and press OK to drain remaining water.
5. Remove the water filter, rinse it thoroughly under running water, and dry it completely.
6. Check the filter assembly direction and reinstall it by placing the center hole over the water-filter guide first.
7. Reinstall the water tank and leave the system clean and dry for storage.

Important

Stop and re-check the setup if the individual has pain, skin pinching, leakage, an unresolved alert, or an insecure connection.