

Carebidet Quick Guide 16

Cleaning & Routine Maintenance

Parts Cleaning and Disinfection

Use this with the training video

This quick guide summarizes the procedure shown in the approved Carebidet training video. For complete operating instructions and cautions, use the Carebidet User Manual and the full video.



Quick steps

1. Disassemble the gender module and detachable cushion from the Diaper Cup.
2. Remove visible soil from the inside and outside of the Diaper Cup using soaked gauze, absorbent cotton, a sponge, or a soft cloth.
3. Apply the approved disinfectant evenly to the inside and outside surfaces and leave it in place for the contact time shown in the training video (10 minutes).
4. Wipe/rinse the disinfected surfaces twice using gauze or a cloth soaked with clean water or sterile saline solution.
5. Repeat the same clean-disinfect-rinse process for the gender module and other approved removable parts.
6. Dry all parts completely, then store them away from direct sunlight in a well-ventilated area.

Important

Stop and re-check the setup if the individual has pain, skin pinching, leakage, an unresolved alert, or an insecure connection.