


Carebidet Quick Guide 18

Fitting & Daily Use

Perineal Care

Use this with the training video

This quick guide summarizes the procedure shown in the approved Carebidet training video. For complete operating instructions and cautions, use the Carebidet User Manual and the full video.



Quick steps

1. Before fitting the cup for perineal care, adjust the Auto Urine Waiting Time from 60 seconds to 20 seconds as shown, then confirm/save the setting.
2. Prepare Carebidet and close the main-body cover before positioning the individual.
3. Fit the Diaper Cup carefully and check that skin is not pinched or caught against the cup or pad.
4. Use the Carebidet controls and the perineal-care procedure shown in the training video for the required care cycle.
5. After care is complete, remove the Diaper Cup carefully and support it on the holder.
6. Re-dress the individual, return the legs/body to a comfortable position, and re-check the skin.
7. Restore any temporary Carebidet settings to the facility/home standard before the next normal use.

Important

Stop and re-check the setup if the individual has pain, skin pinching, leakage, an unresolved alert, or an insecure connection.